

Intended use:

CUANTUM Cura rápido *plus*⁺ is a tissue adhesive intended for wound closure and protection in animals.

Product description:

CUANTUM Cura rápido *plus*⁺ is a high viscosity violet liquid, with a clean and accurate application. In contact with skin and exudates, it solidifies in a few seconds, creating a flexible and resistant layer, sealing the wound and protecting against infection and external agents. Formed transparent film is 100% water resistant and long lasting, being eliminated naturally as affected dermal layer regenerates and grows.

Use instructions:**Prior to application:**

- 1) Clip hair from site if necessary
- 2) Wash with anti-bacterial soap, rinse and dry
- 3) Do not apply any cream or ointment before the adhesive application.

Application:

- 1) Use the product at room temperature.
- 2) Remove the cap unscrewing.
- 3) Approximate the wound edges (if necessary) and apply a thin layer of product on desired area.
- 4) Press container to facilitate adhesive pouring and apply topically on affected area, creating a thin film.
- 5) When applied layer is completely dry to release wound or incision edges (drying time is around 30 seconds – 1 minute, it depends on amount dosis).
- 6) After application, clean the applicator tip with a tissue and screw the cap until the following use.

Removal, storage and conservation:

- Product will slough off naturally as the dermal layer continues to grow.
- Keep at room conditions in a dry place, with device preferably in upright position.
- Adhesive is a fluid liquid. Discard if it has lost fluency.

Precautions and warnings:

- Apply a thin layer. Thick layers delay drying time and may break off prematurely.
- Only external use. Do not apply the product below the epidermis.
- Do not use cellulose, cotton or polyurethane substrates (gauze, cotton-tipped or foam applicators, etc.), as well as basic substances (soaps, etc.). These materials can cause an excessive exothermic reaction.
- Avoid adhesive contact with gloves and other surgical instruments (polymerized adhesive on instruments can be easily removed with acetone). In case of its application in undesired area it can be removed with cellulose.
- Do not use around eyes, burns or infected zones.
- In case of contact with eyes: rinse carefully with water for several minutes. Remove contact lenses, if it is possible.
- Contains cyanoacrylate. It adheres to skin and ayes in a few seconds.
- Keep away from children.
- Only veterinary use



Manufactured by:

CUANTUM MEDICAL COSMETICS, S.L.
Edifici Eureka, Campus UAB, Avda. Can Domènech, s/n
08193 Cerdanyola del Vallès (Barcelona), Spain +34 658 875 835

Rev. 11/2018